

# IMPLANT VARIETY AT A GLANCE

## One-piece implant types

 = length in mm

Abbreviations: m = monotype,  
t = tapered, c = cemented

### Z5m

The trusted implant  
for a variety of indications



### Z5m(t)

The implant for  
immediate implantation  
and soft bone qualities



## Two-piece implant

### Z5c

The successful implant  
for cemented tissue-level  
restoration



	Ø 3,6 mm	Ø 4,0 mm	Ø 5,0 mm
8		Z5m-40-08	Z5m-50-08
10	Z5m-36-10	Z5m-40-10	Z5m-50-10
12	Z5m-36-12	Z5m-40-12	Z5m-50-12

	Ø 4,0 mm	Ø 5,0 mm
8	Z5m(t)-40-08	Z5m(t)-50-08
10	Z5m(t)-40-10	Z5m(t)-50-10
12	Z5m(t)-40-12	Z5m(t)-50-12

	Ø 4,0 mm	Ø 5,0 mm
8	Z5c-40-08	Z5c-50-08
10	Z5c-40-10	Z5c-50-10
12	Z5c-40-12	Z5c-50-12

## Abutments

### Z5c-A00

Straight abutment (5 mm height)



### Z5c-A00-4mm

Straight abutment (4 mm height)

### Z5c-A15

Angled abutment, 15° (5 mm height)



### Z5c-A15-4mm

Angled abutment, 15° (4 mm height)



swiss made 

#### Headquarters

Z-SYSTEMS AG | Werkhofstrasse 5 | CH-4702 Oensingen | Tel. +41 62 388 69 69 | support@zsystems.com

#### US Distribution

Z-SYSTEMS USA Inc. | 284 Monponsett Street, Suite 209 | Halifax, MA 02338 | phone +1 (781) 754-6566 | fax +1 (781) 754-6565 | usa@zsystems.com

www.zsystems.com